

TECH SERVICES DOCUMENT

Shelf Life on ProVent® Fabric

Rev. (11/15)

Shelf life is typically quoted for products or materials which contain components that are expected to degrade and/or lose their performance effectiveness after a period of time.

Kappler has conducted natural aging studies during package validations on surgical gowns made of ProVent® fabric. Testing has been conducted on products that have been through storage and handling and are at least five years old. The testing included the ASTM 1671, Blood and Viral Resistance, and verification of package sterility. Based on these studies, ProVent® retains its performance characteristics for at least five years.

Warning: This information is based on technical data that Kappler believes to be reliable. It is subject to revision as additional knowledge and experience are gained. The website will contain Kappler's most up to date product information, and customers who receive pamphlets, brochures or other literature should be aware that such "hard copy" information may not be as current as the information on Kappler's website. Customers should also recognize that there are uses, environments and chemicals for which Kappler products, garments and/or fabrics are unsuitable. It is the responsibility of the user to review available data and verify that standards. Also, the customer should review all available information on the website to understand the uses- and limitations- on ALL products, garments and fabrics which Kappler makes available.

CAUTION: Do not use for fire protection. Avoid open flame or intense heat.

Kappler Inc, • Technical Services Department • Post Office Box 490 (Shipping: 55 Grimes Dr.) Guntersville, AL 35976 • Toll Free 1-800-600-4019 • 1-256-505-4005 • Fax 1-256-505-4151 • Email: customerservice@kappler.com